

CRIZAL® AVAILABILITY CHART

FOR EYECARE PROFESSIONALS



EFFECTIVE — JUL 2010

Monster scratches appear
on ordinary lenses

Not on Crizal Avancé lenses
with Scotchgard Protector

Crizal Avancé™ lenses
with Scotchgard™ Protector

The most significant breakthrough in
scratch resistance

For more information
www.crizalusa.com



CRIZAL™ AVAILABILITY CHART

for Eyecare Professionals



EFFECTIVE JUL 2010 2010 : Release.04		PLASTIC 1.50			TRIVEX®		Airwear®				THIN&LITE® 1.60, 1.67, 1.74				
		Clear	Xperio® Polarized	Transitions® VI	Trivex®	Trivex® Transitions® VI	Clear	Xperio® Polarized	Transitions VI	Transitions XTRActive	1.60 Clear	1.60 Transitions VI	1.67 Clear	1.67 Transitions VI	Thin&Lite 1.74 ¹
TINTABLE		✓								✓		✓ ²			
PROGRESSIVE LENSES															
NEW	Varilux Physio Enhanced™	•	•	•	•	• G	•	•	•	• G	•	•	•	•	•
	Varilux Physio 360 ^o ™ ³	•	•	•			•	•	•		•	•	•	•	•
	Varilux Physio Short 360 ^o ™ ³						•		•			•	•		
NEW	Varilux® Physio®	•	•	•	•	• G	•	•	•	• G	•	•	•	•	•
NEW	Varilux Physio Short™	•		•	•	• G	•		•		•	•	•	•	
NEW	New Varilux Comfort Enhanced™	•		•			•		•						
NEW	New Varilux Comfort®	•		•			•		•						
NEW	New Varilux Comfort Short™	•		•			•		•						
	Varilux Comfort 360 ^o ™ ³		•					•			•	•	•	• G	
	Varilux Comfort®		•					•			•	•	•	• G	
	Varilux® Panamic®	•	•	•			•	•	•		•	•	•	• G	•
	Varilux Ellipse 360 ^o ™ ³	•		• G			•	•	•		•	• G	•	• G	•
	Varilux® Ellipse®	•		• G			•	•	•		•	• G	•	• G	
DEFINITY	DEFINITY®	•	•	•	•	•	•	•	•		•		•	•	•
	DEFINITY SHORT™	•	•	•	•	•	•	•	•		•		•	•	•
NEW	Ideal™	•	•	•			•	•	•		•	•	•	•	
NEW	Ideal Short™	•	•	•			•	•	•		•	•	•	•	
NEW	Ideal Advanced™	•	•	•			•	•	•		•	•	•	•	
NEW	Ideal Advanced Wrap™	•	•	•			•	•	•			•	•	•	
	Accolade®	•		•	•	•	•	•	•	• G	•	•	•	•	•
	Accolade Freedom™	•		•	•	•	•	•	•	• G	•	•	•	•	•
	SmallFit®	•		• G			•	•	•		•		•	•	
	Ovation®	•		•			•	•	•			•	•		
	Essilor Natural®	•		• G			•		• G		•				
	Adaptar®	•		• G			•								

1 All lenses with Thin&Lite 1.74 material will systematically include choice of New & Improved Crizal® Alizé® or Crizal Avancé™ with Scotchgard™ Protector.
 2 All 1.67 materials can only be tinted in Gray and Brown.
 3 These lenses systematically include a Crizal anti-reflective of your choice.
 G Denotes that this Transitions VI product is available only in Gray.

For more information, please visit www.crizalusa.com

CRIZAL™ AVAILABILITY CHART

for Eyecare Professionals



EFFECTIVE JUL 2010 2010 : Release.04	PLASTIC 1.50			TRIVEX®		Airwear®				THIN&LITE® 1.60, 1.67, 1.74				
	Clear	Xperio® Polarized	Transitions® VI	Trivex®	Trivex® Transitions® VI	Clear	Xperio® Polarized	Transitions VI	Transitions XTRActive	1.60 Clear	1.60 Transitions VI	1.67 Clear	1.67 Transitions VI	Thin&Lite 1.74 ¹
TINTABLE	✓									✓		✓ ²		
SINGLE VISION														
COMING SOON	SV	•	•	•	•	•	•	• G		•	•	•	• G	•
COMING SOON	SV AS				•	•	•	• G		•	•	•	•	•
	Essilor 360 ^{OTM}					•		•				•	• G	•
NEW	Essilor 360° Wrap TM	•	•	•		•	•	•				•	•	
	Essilor Sport Solutions ^{TM 3}						•							
BIFOCALS/ TRIFOCALS														
COMING SOON	FT28	•	•	•	•	•								
	FT28 AS											•		
	FT28 Essilor 360 ^{OTM}	•		•		•								
	FT35	•		•										
	FT45	•												
	ROUND 22	•												
	7x28 (TriFocal)	•		•										
	8x35 (TriFocal)	•		•										
MID DISTANCE / VISUAL FATIGUE LENSES														
	Essilor Anti-Fatigue TM	•				•						•		
	Essilor Computer TM					•								
	Interview TM	•												
	PRIO® Computer Lens	•												

1 All lenses with Thin&Lite 1.74 material will systematically include choice of New & Improved Crizal® Alizé® or Crizal Avancé™ with Scotchgard™ Protector.
 2 All 1.67 materials can only be tinted in Gray and Brown.
 3 Essilor Sport Solutions has six different lens colors available, there is a lens to complement the unique lighting condition of most outdoor sporting activities.
 G Denotes that this Transitions VI product is available only in Gray.

For more information, please visit www.crizalusa.com

CRIZAL™ AVAILABILITY CHART

for Eyecare Professionals



EFFECTIVE JUL 2010 2010 : Release.04		PLASTIC 1.50				Trivex®	MID	POLY		HI-INDEX							
		Clear	Xperio® Polarized	Transitions® VI	Melanin		1.56 Polarized	Clear	Xperio® Polarized	1.60 Clear	1.60 Transitions VI	1.60 Xperio® Polarized	1.67 Clear	1.67 Transitions VI	1.67 Xperio® Polarized		
TINTABLE		✓															
ESSILOR COMPANIES																	
KBco	SV		•						•								•
	SV AS																•
	FT28 Spherical		•						•								
	7x28		•														
	EOST™ PAL		•						•								•
	Fusion I PAL		•														
	Fusion II PAL								•								
Specialty Lens Corporation	Specialty Select™ short corridor					•											
	Specialty Select™ progressive					•											
	iRx Polaroid® SV		•		•	•	•		•		•						•
	iRx Polaroid® SV AS																•
	iRx Polaroid® SV-XL 85mm		•														
	iRx Polaroid® FT28		•		•		•										
	iRx Polaroid® FT35		•		•		•		•								
	iRx Polaroid® 7x28		•		•		•										
	iRx Polaroid® 8x35		•		•		•										
	iRx Polaroid® SHORT				•		•		•								
	iRx Polaroid® PRO		•		•		•				•						•
	iRx Polaroid® SV / PRO Transitions® SOLFX™		•														
Signet Armorlite	Kodak®	•															
	Kodak® Concise™	•				•			•				•				
	Kodak® Precise™	•				•			•				•				
	Kodak® Precise™ Short	•							•				•				

For more information, please visit www.crizalusa.com

CRIZAL™ AVAILABILITY CHART

for Eyecare Professionals



EFFECTIVE JUL 2010 2010 : Release.04		PLASTIC 1.50			Trivex®	POLY			HI-INDEX					Younger Drivewear™ 1.50 / Polycarbonate
		Clear	Polarized	Transitions® VI		Clear	Transitions VI	Polarized	1.60 Clear	1.60 Transitions VI	1.67 Clear	1.67 Transitions VI	1.67 Polarized	
TINTABLE		✓												
NON-ESSILOR COMPANIES														
Shamir	Genesis™	•		•	•				•	•	•	•		
	Genesis w/ Attitude™								•	•				
	Office™	•				•								
	Creation™	•		• G		•			•	•	•	•		
	Piccolo®	•		• G		•			•	•	•	•		
Sola	Compact Ultra™	•		•							•	•		
	VIP®	•		•										
	SOLA MAX®	•		• G										
	Sola One®	•		•							•	•		
AO	Compact™	•		• G				•						
Younger	SV Spherical				•									•
	SV Aspheric				•									
	Image®	•	•	•	•									•
X-Cel	SV			•										
Optima	Optima Resolution™					•								
	Optima Response™					•								
Zeiss	GT2	•		•			•				•	•		
	GT2 Short	•		•			•				•	•		
	Gradal Top	•		•										
Seiko	SFSV										•			
	Proceed II										•			
	Proceed III										•			
Pixel	atLast Composite										•	•		

SPECIALTY ITEMS: on UNCOATED LENSES ONLY

Slaboff, High Adds, Bi-Concave, etc...	•													
--	---	--	--	--	--	--	--	--	--	--	--	--	--	--

G Denotes that this Transitions VI product is available only in Gray.

For more information, please visit www.crizalusa.com



FINISHED SINGLE VISION
Crizal Avancé™ Lenses with Scotchgard™ Protector

Material	Index	Power Range	Size
Thin&Lite® 1.74 (Aspheric)	1.74	-3.00 to -12.00 up to -2.00 cyl	70/75/80mm
		-2.00 cyl available to a -10.00	
Thin&Lite® 1.67 (Aspheric)	1.67	+6.00 to -10.00 up to -2.00 cyl	65/70/75/80mm
		-2.00 cyl available to a -10.00*	
Thin&Lite® 1.67 Transitions® VI Gray (Aspheric)	1.67	-2.00 to -8.00 up to -2.00 cyl	70/75/80mm
		-2.00 cyl available to a -8.00	
Thin&Lite® 1.60 (Aspheric)	1.60	+6.00 to -10.00 up to -2.00 cyl	65/70/75mm
		-2.00 cyl available to a -10.00	
Thin&Lite® 1.60 Transitions® VI Gray (Aspheric)	1.60	+4.00 to -6.00 up to -2.00 cyl	65/70/75mm
		-2.00 cyl available to a -6.00	
Airwear® (Aspheric) & (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70mm
		-2.00 cyl available to a -6.00	
Airwear® Transitions® VI Gray (Aspheric) & (Spherical)	1.59	+4.00 to -6.00 up to -2.00 cyl	65/70mm
		-2.00 cyl available to a -4.00	
Orma® Flat	1.50	+4.00 to -6.00 up to -2.00 cyl	70mm
		-2.00 cyl available to a -4.00	
1.50 Transitions® VI Gray	1.50	+4.00 to -4.00 up to -2.00 cyl	65/70mm
		-2.00 cyl available to a -2.00	

New & Improved Crizal® Alizé®

Material	Index	Power Range	Size
Thin&Lite® 1.67 (Spherical)	1.67	-2.00 to -8.00 up to -2.00 cyl	70/75mm
		-2.00 cyl available to a -6.00	
Thin&Lite® 1.67 (Aspheric)	1.67	+6.00 to -12.00 up to -2.00 cyl	65/70/75/80mm
		-2.00 cyl available to a -10.00*	
Thin&Lite® 1.67 Transitions® VI Gray (Aspheric)	1.67	-2.00 to -8.00 up to -2.00 cyl	70/75/80mm
		-2.00 cyl available to a -8.00	
Thin&Lite® 1.60 (Spherical)	1.60	+6.00 to -10.00 up to -2.00 cyl	65/70/75mm
		-2.00 cyl available to a -10.00	
Airwear® (Aspheric) & (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70mm
		-2.00 cyl available to a -6.00	
Airwear® Transitions® VI Gray (Aspheric) & (Spherical)	1.59	+4.00 to -6.00 up to -2.00 cyl	65/70mm
		-2.00 cyl available to a -4.00	
Orma® (Flat)	1.50	+4.00 to -6.00 up to -2.00 cyl	70mm
		-2.00 cyl available to a -4.00	
1.50 Transitions® VI Gray (Aspheric)	1.50	+4.00 to -4.00 up to -2.00 cyl	65/70mm
		-2.00 cyl available to a -2.00	

*Also available from -10.00 to -12.00 in Spherical only

For more information, please visit www.crizalusa.com

FINISHED SINGLE VISION

Crizal®			
Material	Index	Power Range	Size
Thin&Lite® 1.67 (Spherical)	1.67	-2.00 to -8.00 up to -2.00 cyl	70/75mm
		-2.00 cyl available to a -6.00	
Thin&Lite® 1.67 (Aspheric)	1.67	-2.00 to -8.00 up to -2.00 cyl	75/80mm
		-2.00 cyl available to a -6.00	
Airwear® (Aspheric)	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70mm
		-2.00 cyl available to a -6.00	
Airwear®	1.59	+4.00 to -8.00 up to -2.00 cyl	65/70mm
		-2.00 cyl available to a -6.00	
Orma® Thin	1.50	+4.00 to -5.00 up to -2.00 cyl	70mm
		-2.00 cyl available to a -3.00	
Orma®	1.50	+4.00 to -6.00 up to -2.00 cyl	70/75mm
		-2.00 cyl available to a -4.00	

For more information, please visit www.crizalusa.com